


Connection ISPyB and Scipion: Current status at the ESRF and future plans

Olof Svensson
Data Analysis Unit
ISDD / ESRF



<http://www.esrf.eu/home/news/general/content-news/general/inauguration-of-a-cryo-electron-microscope-platform-at-the-esrf.html>

Computing resources:

- One dedicated server
- 28 cores
- 2 GPU (GTX 1080 Ti)

The screenshot displays the Scipion software interface. The top bar shows the project name 'eput-grid2-launch-on13122017 (opcm01 on cmproc1)' and the user 'ESRF_Debian9 (2017-12-11) d94e7eded1'. The left sidebar contains a tree view of protocols, with 'Import particles' selected. The main area shows a workflow diagram with steps: 'PROJECT', 'scipion - import movies running (done 0/1)', 'motioncorr - motioncorr alignment running (done 1123/1686)', 'gctf - CTF estimation on GPU running (done 1122/1123)', and 'ispyb - monitor to ISPyB at the ESRF running (done 0/1)'. The bottom right window shows the 'run.stdout' log, which contains information about importing files from the ESRF data storage.

```
graph TD
    PROJECT[PROJECT] --> scipion[scipion - import movies running (done 0/1)]
    scipion --> motioncorr[motioncorr - motioncorr alignment running (done 1123/1686)]
    motioncorr --> gctf[gctf - CTF estimation on GPU running (done 1122/1123)]
    gctf --> ispyb[ispyb - monitor to ISPyB at the ESRF running (done 0/1)]
```

run.stdout

```
Square_11423545/Data/FoilHole_13440807_Data_11451488_11451489_2017
1214_1255-14688.mrc
01700: 2017-12-14 13:04:12,026 INFO: Importing file:
/data/visitor/mx2019/cm01/20171211/RAW_DATA/eput-grid2/Images-Disc1/Grid
Square_11423545/Data/FoilHole_13440811_Data_11451488_11451489_2017
1214_1257-14689.mrc
01701: 2017-12-14 13:07:31,308 INFO: Importing file:
/data/visitor/mx2019/cm01/20171211/RAW_DATA/eput-grid2/Images-Disc1/Grid
Square_11423545/Data/FoilHole_13440813_Data_11451488_11451489_2017
1214_1300-14690.mrc
01702: 2017-12-14 13:11:27,380 INFO: Importing file:
/data/visitor/mx2019/cm01/20171211/RAW_DATA/eput-grid2/Images-Disc1/Grid
Square_11423545/Data/FoilHole_13440840_Data_11451488_11451489_2017
1214_1304-14691.mrc
```

Customized scipion startup script:

- Creates a new Scipion project and executes it.
- Command line driven, started by the CM01 scientists (GUI to be developed).
- Uses Scipion un-modified actors for movie import, motioncorr2 and Gctf.
- ISPyB monitor developed using the one written by (Kevin Savage and Jose Miguel de la Rosa Trevin) as base.
- Runs very reliable! Very nice features like possibility of individual start and stop of actors.
- Uploads meta-data to ISPyB using three web service calls.
- In operation since November 2017, 7 user experiments so far, > 25000 movies uploaded to ISPyB.

<https://github.com/antolinos/ispyb-client/wiki/EM-Workflow>

addMovie :

- proposal: Proposal type and proposal number of the current experiment. For instance: mx415
- sampleAcronym: Acronym of the sample.
- movieDirectory: [DATACOLLECTION.imageDirectory] Directory that contains movies and will identify the grid square as a data collection.
- movieFullPath: [MOVIE.movieFullPath] The full path of the movie.
- movieNumber: [MOVIE.movieNumber] The unique number of a movie for a grid square.
- micrographFullPath: [MOVIE.micrographFullPath] Full path of a micrograph produced by instrument.
- micrographSnapshotFullPath: [MOVIE.micrographSnapshotFullPath] Full path of a snapshot of the micrograph.
- xmlMetaDataFullPath: [MOVIE.xmlMetaDataFullPath] Full path of a XML file containing movie metadata.
- voltage: [DATACOLLECTION.voltage] Microscope voltage in kV.
- sphericalAberration: [DATACOLLECTION.sphericalAberration] Spherical aberration in mm.
- amplitudeContrast: [DATACOLLECTION.amplitudeContrast] Amplitude contrast.
- magnification: [DATACOLLECTION.magnification] Magnification.
- pixelSizeOnImage: [DATACOLLECTION.pixelSizeOnImage] scanned pixel size in microns.
- numberOfImages: [DATACOLLECTION.numberOfImages] Number of frames in a movie.
- dosePerImage: [MOVIE.dosePerImage] Dose per image in $e/\text{\AA}^2$.
- positionX: [MOVIE.positionX] X position of movie.
- positionY: [MOVIE.positionY] Y position of movie.
- beamlineName
- startTime SimpleDateFormat dt = new SimpleDateFormat("yyyy-mm-dd hh:mm:ss");

addMotionCorrection

- proposal: Proposal type and proposal number of the current experiment. For instance: mx415
- movieDirectory: Directory that contains movies and will identify the grid square as a data collection.
- micrographFullPath: Full path of a MRC file containing the micrograph.
- correctedDoseMicrographFullPath: Full path of a MRC file containing the micrograph corrected for dose.
- micrographSnapshotFullPath: Full path of a PNG file containing snapshot of micrograph.
- driftPlotFullPath: Full path of a PNG file containing plots of global shifts.
- logFileFullPath: Full path to the log file of motion correction.

addCTF

- proposal: Proposal type and proposal number of the current experiment. For instance: mx415
- movieFullPath: [MOVIE.movieFullPath] The full path of the movie.
- spectralImageSnapshotFullPath: Full path of a JPEG snapshot of the spectra image.
- spectralImageFullPath: Full path to the MRC file of the spectra image.
- defocusU: Defocus U in Å.
- defocusV: Defocus V in Å.
- angle: angle in degrees.
- crossCorrelationCoefficient: The Cross Correlation Coefficient (CCC).
- resolutionLimit: Resolution limit in Å.
- estimatedBfactor: Estimated B factor in Å².
- logFilePath: full path to the log file of CTF estimation

PRESENTATION OF RESULTS IN EXI

EH 23-11-2017 16:37:53

/data/visitor/mx415/cm01/20171123/RAW_DATA/testsecretin-grid1/Images-Disc1/GridSquare_23722826/Data

Sample	Sample6
Grid Squares	1
Voltage	300000
Spherical Aberration	2.7
# Frames	40
Amplitude Contrast	10%
Scanned Pixel Size	



2017-11-23 16:37:52

Movies

7

Motion Correction

100.0%

CTF

100.0%

Voltage 300000

Magnification 130000

Directory /data/visitor/mx415/cm01/20171123/RAW_DATA/testsecretin-grid1/Images-Disc1/GridSquare_23722826/Data

Movie #3385 30-11-2017 10:55:23

Movie


Motion Correction

CTF

Number 3385

X, Y -0.000, 0.000

Dose 0.1300





Total Motion 11.3

Avg. Motion/frame 0.3

Frame Range 1 - 40

Dose 1.0000



Resolution Limit 2.902




Correlation 0.068865

Defocus U 22280.72

Defocus V 21719.82

Angle 62.99

Estimated B factor 86.87




Extend automatic processing

- Particle picking
- Particle extract
- 2D classification

Mail-in operation






- Not possible to give users access to Cryo-EM data acquisition software like for MXCuBE
- Users will choose remotely which grid squares to use for data collection



DATA MODEL



11



Page 11 | Connecting ISPyB and Scipion | Olof Svensson | 2018-01-30

5



ACKNOWLEDGEMENTS

- Isai Kandiah, Alejandro de Maria, Solange Delagenière, Christoph Mueller-Dieckmann and Gordon Leonard (ESRF)
- Michael Hons (EMBL Grenoble)
- Gregory Effantin (IBS)
- Kevin Savage, Karl Levik, Neil Smith, Stuart Fischer and Alun Ashton (DLS)
- Pablo Conesa, Jose Miguel de la Rosa Trevin and Yaiza Rancel Gil (CNB)
- The ESRF Data Analysis Unit, Beamline Control Unit and Structural Biology group